

Culture broths for blood and cerebrospinal fluid

CULTURE BROTHS FOR BLOOD AND CEREBROSPINAL FLUID**Ref.No.****Packaging****Volume**

Culture broths for blood and cerebrospinal fluid allow isolation of microorganisms and secure the drawn samples for further diagnostics procedures.

Hemomedium

For the enrichment of aerobic bacteria from blood.

5001

vial

50 ml

5002

vial

15 ml

5003

vial

5 ml

Anaeromedium

For the enrichment of anaerobic bacteria from blood.

5004

vial

50 ml

5005

vial

15 ml

5006

vial

5 ml

Meningomedium

For the enrichment of microorganisms from cerebrospinal fluid.

5009

vial

50 ml

5008

vial

15 ml

5007

vial

5 ml