

Harmonised Method EP/USP/JP

CULTURE MEDIA ACCORDING TO HARMONISED METHOD EP/USP/JP	Ref.No.	Packaging	Volume
Cetrymide Agar For the selective isolation of <i>Pseudomonas aeruginosa</i> .	9001 50	powder	500 g
	9001 05	bottle	500 ml
	9001 02	bottle	200 ml
	9001 01	bottle	100 ml
	9001 04	plate	90 mm
Columbia Agar For the cultivation of anaerobic microorganisms.	9002 50	powder	500 g
	9002 05	bottle	500 ml
	9002 02	bottle	200 ml
	9002 01	bottle	100 ml
	9002 04	plate	90 mm
MacConkey Agar with crystal violet For the isolation of <i>Escherichia coli</i> .	9003 50	powder	500 g
	9003 05	bottle	500 ml
	9003 02	bottle	200 ml
	9003 01	bottle	100 ml
	9003 04	plate	90 mm
MacConkey Broth For the selective enrichment of <i>Escherichia coli</i> .	9004 50	powder	500 g
	9004 05	bottle	500 ml
	9004 02	bottle	200 ml
	9004 01	bottle	100 ml
	9004 11	tube	10 ml
Mannitol Salt Agar (MSA) For the isolation of <i>Staphylococcus aureus</i> .	9005 50	powder	500 g
	9005 05	bottle	500 ml
	9005 02	bottle	200 ml
	9005 01	bottle	100 ml
	9005 04	plate	90 mm
Mossel Broth For the selective enrichment of <i>Enterobacteriaceae</i> .	9006 50	powder	500 g
	9006 05	bottle	500 ml
	9006 02	bottle	200 ml
	9006 01	bottle	100 ml
	9006 11	tube	10 ml
Potato Dextrose Agar For the isolation of yeasts and molds.	9007 50	powder	500 g
	9007 05	bottle	500 ml

	9007 02	bottle	200 ml
	9007 01	bottle	100 ml
	9007 04	plate	90 mm
<hr/>			
R2 Agar	9018 50	powder	500 g
For the isolation and enumeration of bacteria from water.	9018 05	bottle	500 ml
	9018 02	bottle	200 ml
	9018 01	bottle	100 ml
	9018 04	plate	90 mm
<hr/>			
Reinforced Clostridial Medium Agar (RCM)	9008 50	powder	500 g
For the isolation of anaerobic microorganisms.	9008 05	bottle	500 ml
	9008 02	bottle	200 ml
	9008 01	bottle	100 ml
	9008 04	plate	90 mm
<hr/>			
Rappaort - Vassiliadis (RSV) Broth	9009 50	powder	500 g
For the selective enrichment of <i>Salmonella</i> spp.	9009 05	bottle	500 ml
	9009 02	bottle	200 ml
	9009 01	bottle	100 ml
	9009 11	tube	10 ml
<hr/>			
Sabouraud Dextrose Agar	9010 50	powder	500 g
For the isolation and enumeration of fungi.	9010 05	bottle	500 ml
	9010 02	bottle	200 ml
	9010 01	bottle	100 ml
	9010 04	plate	90 mm
<hr/>			
Sabouraud Dextrose Broth	9011 50	powder	500 g
For the enrichment of yeasts and molds.	9011 05	bottle	500 ml
	9011 02	bottle	200 ml
	9011 01	bottle	100 ml
	9011 11	tube	10 ml
<hr/>			
Tryptic Soy Agar (TSA)	9012 50	powder	500 g
For the cultivation of a wide variety of microorganisms.	9012 05	bottle	500 ml
	9012 02	bottle	200 ml
	9012 01	bottle	100 ml
	9012 04	plate	90 mm
<hr/>			
Tryptic Soy Broth (TSB)	9013 50	powder	500 g
For the enrichment of a wide variety of microorganisms.	9013 05	bottle	500 ml
	9013 02	bottle	200 ml
	9013 01	bottle	100 ml
<hr/>			
Violet Red Bile Glucose (VRBG) Agar	9014 50	powder	500 g
For the isolation of <i>Enterobacteriaceae</i> .	9014 05	bottle	500 ml
	9014 02	bottle	200 ml
	9014 01	bottle	100 ml
	9014 04	plate	90 mm

Xylose Lysine Dextrose (XLD) Agar	9016 50	powder	500 g
For the isolation of <i>Salmonella</i> spp. and <i>Shigella</i> spp.	9016 05	bottle	500 ml
	9016 02	bottle	200 ml
	9016 01	bottle	100 ml
	9016 04	plate	90 mm
Buffered Sodium Chloride-Peptone Solution	9017 50	powder	500 g
For the microbiological examination of non-sterile products.	9017 05	bottle	500 ml
	9017 02	bottle	200 ml
	9017 01	bottle	100 ml
Buffered Sodium Chloride-Peptone Solution and Tween 80, 1g/l	9019 05	bottle	500 ml
For the microbiological examination of non-sterile products.	9019 02	bottle	200 ml
	9019 01	bottle	100 ml