

Supplements

SUPPLEMENTS	Ref.No.	Packaging	Volume
Biovitex Growth factor for the cultivation of a wide variety of fastidious microorganisms. <i>Chocolate Agar base (GC Agar) (Ref.No. 109)</i>	420009	vial (lyophilised)	10x500 ml
BCYE (Cysteine, α-ketoglutarate) Growth supplement for the isolation of <i>Legionella spp.</i> <i>BCYE Agar (Ref.No. 8049)</i>	X196	vial (lyophilised)	10x500 ml
Campylobacter supplement (Amphotericin, Cefsulodin, Trimethoprim, Vancomycin, Polymyxin) Selective supplement for the isolation of <i>Campylobacter spp.</i> <i>Campylobacter Agar with 10% sheep blood (Ref.No. 1291)</i>	3055	vial (fluid)	10 ml
Chloramphenicol Selective supplement for the isolation of yeasts and molds. <i>YGC Agar base (Ref.No. 287)</i>	X009	vial (lyophilised)	10x500 ml
CIN supplement (Cefsulodin, Irgasan, Novobicin) Selective supplement for the isolation of <i>Yersinia spp.</i> <i>Yersinia CIN Agar base (Ref.No.123)</i>	X120	vial (lyophilised)	10x500 ml
Cefoperazone, Amphotericin Selective supplement for the isolation of <i>Campylobacter spp.</i> <i>mCCDA base according to ISO 10272-1:2006 (Ref.No. 281)</i>	X112	vial (lyophilised)	10x500 ml
Clostridium supplement (Cycloserine, Cefoxitin) Selective supplement for the isolation of <i>Clostridium spp.</i> <i>Clostridium Difficile Agar (Ref.No.1005)</i>	3057	vial (fluid)	10 ml
CN supplement (Cetrymide, Nalidixic Acid) Selective supplement for the isolation of <i>Pseudomonas aeruginosa</i> . <i>Pseudomonas CN Agar base (Ref.No. 238)</i>	X107	vial (lyophilised)	10x500 ml
Cycloserine Selective supplement for the isolation of <i>Clostridium perfringens</i> . <i>TSC Agar base (Tryptose Sulphite Cycloserine Agar) (Ref.No. 283)</i>	X194	vial (lyophilised)	10x500 ml
Cycloserine, Cefoxitin	X093	vial (lyophilised)	10x500 ml

Selective supplement for the isolation of *Clostridium perfringens* from clinical samples.

C.V.T.C (Cefoperazone, Vancomycin, Trimethoprim, Cycloheximide)

Selective supplement for the isolation of *Campylobacter spp.* from food and environmental samples.

Campylobacter Enrichment Broth base (Bolton Broth) (Ref.No. 289)

X131 vial (lyophilised) 10x500 ml

Egg Yolk Emulsion

A sterile emulsion of egg yolks for use in bacteriological culture media.

Bacillus Cereus Agar base (M.Y.P) (Ref.No. 265)

X073 bottle 100 ml

Egg Yolk Tellurite Emulsion

A sterile egg yolk containing potassium tellurite for use in bacteriological culture media.

Baird Parker Agar base (Ref.No. 131)

X085 bottle 100 ml

Fraser supplement (5% Ferric Ammonium Citrate)

Selective supplement for the isolation of *Listeria spp.*

Fraser Broth base (Ref.No. 223)

3059 vial (fluid) 10 ml

Gentamycin

Selective supplement for the isolation of yeasts and molds.

Sabouraud Dextrose Agar (Ref.No.114)

3063 vial (fluid) 10 ml

HTM (Haemophilus Test Medium) supplement (Factor X and V)

Enrichment supplement for the isolation of *Haemophilus spp.*

HTM (Haemophilus Test Medium) Agar base (Ref.No. 3506)

3127 vial (fluid) 10x500 ml

Kanamycin 75 mg/l

Selective supplement for the isolation of *Clostridium spp.*

Schaedler Agar base (Ref.No.138)

X018 vial (lyophilised) 10x500 ml

Colistin, Nalidixic Acid

Selective supplement for the isolation of gram positive cocci.

Columbia Agar base (Ref.No.107)

X011 vial (lyophilised) 10x500 ml

Nalidixic Acid

Selective supplement for the isolation of *Listeria monocytogenes*.

X072N vial (lyophilised) 10x500 ml

Legionella GVPC

Selective supplement for the isolation of *Legionella spp.*

GVPC Agar (Charcoal Yeast Extract Agar with selective supplements) (Ref.No. 8051)

X195 vial (lyophilised) 10x500 ml

Neomycin 75

Selective supplement for the isolation of *Clostridium spp.*

X015 vial (lyophilised) 10x500 ml

Novobiocin

Selective supplement for the isolation of *Salmonella spp.*

Mueller Kauffman Tetrathionate Broth (MKTtn) (Ref.No. 282)

X150 vial (lyophilised) 10x500 ml

<p>Oxacillin, Polymyxin B Selective supplement for the isolation of methicillin resistant <i>Staphylococcus aureus</i> (MRSA). <i>O.R.S.I.M base (Oxacyllin Resistance Staphylococci Isolation Medium) (Ref.No. 269)</i></p>	X192	vial (lyophilised)	10x500 ml
<p>Oxytetracycline Selective supplement for the isolation of yeasts and molds. <i>OGYE Agar base (Ref.No. 199)</i></p>	X089	vial (lyophilised)	10x500 ml
<p>Oxford Listeria supplement (Cycloheximide, Colistin, Acryflavine, Cefofetan, Fosfomycin) Selective supplement for the isolation of <i>Listeria monocytogenes</i>. <i>Oxford Agar base (Ref.No. 236)</i></p>	X122	vial (lyophilised)	10x500 ml
<p>Palcam supplement (Polymyxin, Acryflavine, Ceftazidime) Selective supplement for the isolation of <i>Listeria monocytogenes</i>. <i>Palcam Agar base (Ref.No. 230)</i></p>	X144	vial (lyophilised)	10x500 ml
<p>Polymyxin B Selective supplement for the isolation of <i>Bacillus cereus</i>. <i>Bacillus cereus Agar base (M.Y.P) (Ref.No. 265)</i></p>	X193	vial (lyophilised)	10x500 ml
<p>PolyVite X Growth factor for the cultivation of a wide variety of fastidious microorganisms. <i>Chocolate Agar base (GC Agar) (Ref.No. 109)</i></p>	3132	vial (fluid)	10 ml
<p>RPF (Bovine fibrinogen, Rabbit plasma, Trypsin inhibitor, Potassium tellurite) Selective supplement for the isolation of <i>Staphylococcus aureus</i>. <i>Baird Parker Agar + RPF (Ref.No. 3226)</i></p>	X086	vial (lyophilised)	10x500 ml
<p>Semi Fraser Supplement (5% Ferric Ammonium Citrate) Selective supplement for the isolation of <i>Listeria spp.</i> <i>Semi-Fraser Broth base (Ref.No. 226)</i></p>	3062	vial (fluid)	10 ml
<p>Ticarcillin, Polymyxin Selective supplement for the isolation of <i>Burkholderia (Pseudomonas) cepacia</i>. <i>Burkholderia Cepacia Agar base (Ref.No. 221)</i></p>	X140	vial (lyophilised)	10x500 ml
<p>Vancomycin, Nalidixic Acid Selective supplement for the isolation of gram negative anaerobes from clinical samples.</p>	X090	vial (lyophilised)	10x500 ml