

# Transport and growth media (DIPSLIDE)

TRANSPORT AND GROWTH MEDIA (DIPSLIDE)	Ref.No.	Packaging
<p><b>MYCAGAR</b> (Fungisel Agar/Sabouraud Dextrose Chloramphenicol Agar) For the isolation of fungi.</p>	<b>7050</b>	10 slides
<p><b>URAGAR</b> (CLED Agar/MacConkey Agar with crystal violet) For the isolation and differentiation of bacteria from urine.</p>	<b>7040</b>	10 slides
<p><b>VRBG</b> Violet Red Bile Glucose Agar / <b>PCA</b> Plate Count Agar For detection of <i>Enterobacteriaceae</i> on surfaces and in liquids.</p>	<b>7031</b>	10 slides
<p><b>Tryptic Soy Agar (TSA) with TTC / Rose Bengal Agar</b> For detection of yeasts and molds on surfaces and in liquids.</p>	<b>7032</b>	10 slides